



CARDIOVASCULAR ASSOCIATES OF AMERICA

Sample Equipment List

Setting up your clean procedure room correctly is the first step in delivering safe and effective vein care. This checklist outlines the core equipment and room requirements needed to support superficial laser, radiofrequency, and vein closure procedures. Practices can use it as a quick-start reference when building or converting a space for vein treatment.

Equipment for Sterile Procedure Room - Vein Clinic

Superficial Laser/Radiofrequency Procedures/Vein Closure Procedures Performed by MDs

Tilt Table or Stretcher

Ultrasound - vein mapping/access

Wall Monitor/TV Monitor

Additional lighting, such as LED surgical-grade lighting

Instrument Tray/Mayo Stand

Hazardous Waste Receptacle

Endovascular Laser/Generator or Radiofrequency Unit/Generator

Injector Pump (tumescent)

Autoclave if not using disposable packs

Small undercounter refrigerator

Room Notes: The stress exam room must be at least equivalent in size to a Stress Echo Room

Must have a dimmer for overhead light

Sink

Doors

Treated and cleaned like a clean room rather than a sterile procedure room

~ Contact Us for More Information ~

matt@strategyco.com • kaitlin.scholz@strategyco.com • brown@cvausa.com